

# NOVABONE<sup>®</sup>

BIOACTIVE SYNTHETIC GRAFT

## NovaBone<sup>®</sup> Morsels

OS-Si<sup>+</sup> Morsels provide users with a bone graft substitute that signals and recruits bone forming cells while providing a matrix of open, interconnected pores for accelerated bone formation to occur.

NovaBone Morsels are designed as a continuous porous structure comparable to cancellous bone. The open structure allows for rapid vascularization and mineralization and permits complete device absorption and replacement by new bone.

### ENGINEERED BONE REGENERATION



OS-Si<sup>+</sup> Morsels are available in a sealed cup or pre-loaded syringe with a Leur Lock attachment which allows for easy expression of blood directly to the morsels for improved handling characteristics during implantation.

#### NovaBone Morsels

Reference No.	Description
NB0805	5cc; 2-5mm, NovaBone Morsels
NB0810	10cc; 2-5mm, NovaBone Morsels
NB0815	15cc; 2-5mm, NovaBone Morsels
NB1805	5cc; 1-2mm, NovaBone Morsels
NB1810	10cc; 1-2mm, NovaBone Morsels
NB1815	15cc; 1-2mm, NovaBone Morsels

#### NovaBone Morsels – Syringe Delivery

Reference No.	Description
NB3810	10cc; 2-5mm, NovaBone Morsels
NB3815	15cc; 2-5mm, NovaBone Morsels
NB4810	10cc; 1-2mm, NovaBone Morsels
NB4815	15cc; 1-2mm, NovaBone Morsels